



E-Prep® BL-110

Burnishing Liquid

BL-110	is a very effective low foaming burnishing liquid concentrate. It is used to obtain a bright clean surface on brass, copper, bronze and steel.
BL-110	is clean rinsing and leaves no tacky residual film. Parts burnished with BL 110 will resist tarnishing for extended periods of time.
BL-110	is applicable to vibratory or rotary barrel type equipment.
BL-110	is biodegradable.

Operating Instructions

Eight (8) fluid ounces (1/2 pint) of **BL-110** concentrate is dissolved in 2 gallons of water. The diluted concentrate is added as required to obtain the desired level of luster.

BL-110 operates best at temperatures from 60°F to 120°F.

IMPORTANT NOTICE! For Industrial Use Only

The following is made in lieu of all warranties, expressed or implied, including the implied warranties of merchantability and fitness for purpose: seller's and manufacturer's only obligation shall be to replace such quantity of the product as proved to be defective. Before using, user shall determine the suitability of the product for its intended use, and user assumes all risk and liability whatsoever in connection therewith. Neither seller nor manufacturer shall be liable either in tort or in contract for any loss or damage, direct, incidental or consequential, arising out of the use or the inability to use the product.